


## Pressure Injury Abatement Solutions




Our consumables and pressure management products help clinical staff to maximize patient safety by reducing shear/friction forces associated with pressure injuries.

# Pressure Equalization Pads

## Tempur-Pedic® For Optimized Patient Care

As many as 2.5 million<sup>1</sup> acute care patients suffer pressure injuries in the US each year, with many occurring in orthopedic and spine surgeries. The total cost of treatment is an astounding \$11 billion annually.<sup>1</sup>

Our exclusive Tempur-Pedic® memory foam pads and positioners provide support, positioning and pressure redistribution during surgical procedures. Studies show Tempur-Pedic material offers up to four times the pressure reduction of conventional polyurethane in controlled testing.<sup>2</sup>



### Features

1. Moisture impenetrable, vapor permeable membrane
2. Flame resistant barrier mattress assembly
3. Tempur 110kg/cubic meter density for peak loading force distribution
4. Tempur-HD soft, dense absorption reduces localized pressure concentration

1. "Hospital Acquired Pressure Injuries HAPU(I) Prevention" Joint Commission Center for Transforming Healthcare

2. Support Surfaces: Viscoelastic Expanded Urethane versus Conventional Foam. Eric Flam, PhD, PE and Loretta Raab, RPh, CCP, ECPM Jan/Feb 2005, pp 29 – 34



## ProneView® Protective Helmet System

The ProneView Protective Helmet System was engineered to optimize pressure redistribution and reduce the risk of central retinal artery occlusion by minimizing pressure on the eyes. The mirror base allows continuous eye monitoring during surgery, so your anesthesiologist can perform all necessary checks with ease. It's also adjustable to accommodate kyphotic patients.

### Features

- A proven system of adjustable stability
- Superior fit, pressure distribution and eye protection
- Quick and easy adjustments
- Superior protection and versatility
- Cushions are available in small, medium and large sizes and vacuum sealed for easy storage.
- More than 4 million cases worldwide




## GentleTouch® Prone Foam Pillows

Keeps the patient's head stable and comfortable during prone procedures. Helps equalize facial pressure while safely supporting the head in the prone position.

### Features

- Single-use GentleTouch Prone Positioning Pillows support the patient's head for safe prone positioning, without hindering access to the anesthesia circuit.
- Made from our exclusive Contourethane™ to equalize pressure on the face and Supportethane™ to cradle the head for maximum safety and comfort.
- There's no glue anywhere on the inside edge, which allows the foam to more easily conform to the patient's face.



## Patient Care Kits Protect Your Patients & Table

Mizuho OSI Patient Care Kits provide support, positioning and isolation during surgical procedures in the prone, lateral or supine position. The vacuum-packed packaging also saves on storage.

### Features

- We optimized our disposable drapes and positioners for use with our specialty tables for maximum performance.
- With over 100 products, you are sure to find just what you need.
- Our extensive product selection includes arm cradles, pad covers, pillows, headrests, perineal post covers, boot liners and more.