

Levó™

Specifications & Components

Specifications

Patient Weight Capacity 650 lbs. (295 kgs)

Adjustment Ranges

Lateral – Dual Motion Base 7.6 in. (19.3 cm)
 Lateral - Lateral Base 6 in. (15.2 cm)
 Longitudinal 6.25 in. (15.9 cm)
 (Dual Motion Base only)
 Height adjustment 6.7 in. (17 cm)
 (Z Tower only)

Arm Ranges

Longitudinal 7.25 in. (18.4 cm)
 Vertical 14.5 in. (36.8 cm)
 Yaw (Ball Joint) ±40°
 Roll (Ball Joint) ±20°

Adjustable Mirror

Lateral tilt ±45°
 Longitudinal tilt -10° to +45°



Ordering Information

Trios® / STS Base Packages

REF 7887-400	Levó Premium Package
7887-300	Levó Standard Package
7887-200	Levó Cervical Package
7887-100	Levó Prone Package

Modules

REF 7887-010	Levó Dual Motion Base
7887-020	Levó Lateral Base
7887-050	Levó Arm
7887-080	Levó Control Handles
7887-840	Contoured Chest Pad
7887-831	Raised Chest Support Plate
7887-070	Levó Prone Platform
7887-060	Levó Z Tower
7887-801	QuickConnect Aluminum Skull Clamp
7887-802	QuickConnect Radiolucent DORO®
7887-803	QuickConnect Radiolucent Mayfield®

Storage Cart

REF 7887-090	Levó Storage Cart
---------------------	-------------------

Levó is ONLY approved and validated for use with the Mizuho OSI Trios® Table System and STS. Use of the device with any other system may result in harm to the patient, the table top, the device or the healthcare professional.

Package Module Components

Module	Prone	Cervical	Standard	Premium
Arm	•	•	•	•
Dual Motion Base		•	•	•
Lateral Base	•			
QuickConnect Aluminum Adapter		•	•	•
Prone Platform	•		•	•
Control Handles		•	•	•
Z-Tower				•
Contoured Chest Pad		•	•	•
Raised Support Plate		•	•	•



Mizuho OSI reserves the right to make changes in specifications and/or to discontinue any product at any time without notice. Mizuho OSI is a Delaware Corporation. Manufactured in the USA.

