



new

Levó™  
Sterile Drape

Visualization With  
Interoperative Control




# Levó™ Sterile Drape

The Levó Head Positioning System introduces a new level of control and precision to head positioning for all spine procedures from cervical to sacrum.

The innovative sterile drape works as a system with Levó to allow full visualization of and interoperative access to all Levó controls without breaking the sterile field.

- Considerable time savings in draping and interoperative adjustments
- Quick and precise articulation without going under the drape
- No need for the surgeon to hold the skull clamp while waiting for critical adjustments
- Surgeon controls all adjustments without breaking the sterile field



1. Sterile drape
2. Clear window for visibility of Levó
3. Sterile sleeves for Levó Control Handles
4. Surgical area opening with expanded outline
5. Hook & Loop Tabs (4)

## Levó Protective Arm Sleeve\*

Also available for Levó is a non-sterile Protective Arm Sleeve that protects the Levó Arm from fluids during surgical procedures.



## Ordering Information

Part Number	Description	Bulk	Protective Arm Sleeve	ShearGuard Chest Pad Cover	Hip/Thigh Pad Cover	Foam Arm Cradles	ProneView Cushion	GentleTouch Cushion	Sterile Drape
7893	Cervical Sterile Drape with Clear Window	6/cs							•
7892	Protective Arm Sleeve	12/cs	•						
7890	Non-sterile Cervical, SG	6/cs	•	•	•				
7887PV	Prone Platform Pkit with ProneView®	6/cs	•	•	•	•	•		
7887GT	Prone Platform Pkit with GentleTouch®	6/cs	•	•	•	•		•	



30031 Ahern Avenue  
 Union City, CA 94587-1234 USA  
 Telephone: 510-429-1500  
 Toll Free: 800-777-4674  
 Fax: 510-429-8500  
 Outside USA: +1-510-429-1500  
[mizuho.com](http://mizuho.com)  
[newhipnews.com](http://newhipnews.com)

**Patient Care Kits Include:** Levó Protective Arm Sleeve, ShearGuard® Chest Pad Cover and Hip and Thigh Pad Covers. Levó Cervical Sterile Drape may be used with all Levó Patient Care Kits.

Note: Mizuho OSI is constantly improving its products. All specifications are subject to change without notice. Mizuho OSI is a Delaware Corporation. Manufactured in the USA. Patents Pending.